

Release Notes

Dentsply Sirona Implants Atlantis[®] suprastructures FLO-S library for 3Shape

Version 3.0.0

- Library is extended to support the new EV connection for DS PrimeTaper and AstraTech Implant EV.
- Library is extended to support Elos Accurate[®] Hybrid Base Non-Engaging and Hybrid Base H Non-Engaging. Please refer to compatibility chart for further information.

Version 2.0.0

- Library is extended and includes interfaces for Nobel Biocare and Straumann.

Version 1.0.1

- Correction of the naming for EV 5.4.

Version 1.0.0

- The library is used for digital design guidance for Atlantis suprastructures and for building and printing models.
- Available for 3Shape software version 2018 and 2019 (and higher.)