

## Release Notes

# Dentsply Sirona Implants Atlantis® suprastructures FLO-S library for Exocad

## Version 3.0.0

- Library is extended to support the new EV connection for DS PrimeTaper and AstraTech Implant EV.
- Library is extended to support Elos Accurate® Hybrid Base Non-Engaging and Hybrid Base H Non-Engaging. Please refer to compatibility chart for further information.

## Version 2.0.0

- Library is extended and includes interfaces for Nobel Biocare and Straumann.

## Version 1.0.1

- Correction of the naming for EV 5.4.

## Version 1.0.0

- The library is used for digital design guidance for Atlantis suprastructures and for building and printing models.