

Temp Abutment Profile EV library compatibility chart

This chart describes the content of the Dentsply Sirona Temp Abutment Profile EV library. This library is used when designing restorations on Temp Abutment Profile EV. The library can be used for 3Shape and Exocad design software.

					
Temp Abutment Profile EV	Atlantis FLO (ref)	Atlantis IO FLO (ref)	Elos Accurate® IO Scan Body	Elos Accurate® ScanBody22™	Printed Model Analog - EV / Elos Accurate® Analog for Printed Models
Temp Abutment Profile EV 4.2 (yellow) (26256)	P-06 (34598)	IO-P-06 (35248)	IO 3C-A	SB-ATE42P-1	68020088 / PMA-ATE42P-1
Temp Abutment Profile EV 4.8 (blue) (26257)	P-07 (34599)	IO-P-07 (35249)	IO 3C-B	SB-ATE48P-1	68020089 / PMA-ATE48P-1

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

Note: The product assortment can vary from country to country. Please check availability with your local Customer Service Representative. TempAbutment Profile EV is not available in digital libraries in US.