

Astra Tech Implant EV PrimeTaper EV Implant OmniTaper EV Implant

TitaniumBase EV ASA library compatibility chart

This chart describes the content of the Dentsply Sirona TitaniumBase EV ASA library. This library is used when designing restorations on TitaniumBase EV ASA. The library can be used for 3Shape and Exocad design software.



TitaniumBase EV ASA	ScanBody EV / Elos Accurate® ScanBody22™	Atlantis FLO* (ref)	Atlantis IO FLO* (ref)	Elos Accurate® IO Scan Body	3Shape / CMT Scan Body	Printed Model Analog - EV / Elos Accurate® Analog for Printed Models
TitaniumBase EV ASA (XS) - green 68012550 - 68012553 TitaniumBase EV ASA (XS) NI - green 68012570 - 68012573	68020517 / SB-ATE30-1	P-01 (34593)	IO-P-01 (35243)	IO 3B-A	DEG	68020083 / PMA-ATE30-1
TitaniumBase EV ASA (S) - purple 68012554 - 68012557 TitaniumBase EV ASA (S) NI - purple 68012574 - 68012577	68020518 / SB-ATE36-1	P-02 (34594)	IO-P-02 (35244)	IO 3B-B	DEP	68020084 / PMA-ATE36-1
TitaniumBase EV ASA (M) - yellow 68012558 - 68012561 TitaniumBase EV ASA (M) NI - yellow 68012578 - 68012581	68020519 / SB-ATE42-1	P-03 (34595)	IO-P-03 (35245)	IO 3B-C	DEY	68020085 / PMA-ATE42-1
TitaniumBase EV ASA (L) - blue 68012562 - 68012565 TitaniumBase EV ASA (L) NI - blue 68012582 - 68012585	68020520 / SB-ATE48-1	P-04 (34596)	IO-P-04 (35246)	IO 3B-D	DEB	68020086 / PMA-ATE48-1
TitaniumBase EV ASA (XL) - brown 68012566 - 68012569 TitaniumBase EV ASA (XL) NI - brown 68012586 - 68012589	68020521 / SB-ATE54-1	P-05 (34597)	IO-P-05 (35247)	IO 3B-E	DE5	68020087 / PMA-ATE54-1

* Atlantis FLO and Atlantis IO FLO can only be used with indexed TitaniumBase EV ASA. Use Atlantis FLO-S with non-indexed (NI) components instead:
TitaniumBase EV ASA (XS) NI: A04K (68020018) TitaniumBase EV ASA (S) NI: A04L (68020019) TitaniumBase EV ASA (M) NI: A04M (68020020)
TitaniumBase EV ASA (L) NI: A04N (68020021) TitaniumBase EV ASA (XL) NI: A04O (68020022)

All trademarks and company names are the property of their respective owners. Elos Accurate® is a registered trademark.

Note: Products are available from the respective manufacturers and their distributors. The product assortment can vary from country to country. Please check availability with your local Customer Service Representative.

