

Release notes DE-2504

Dentsply Sirona Pre-milled Blanks Library for Exocad

Version 1.0.0

- The library is used in the construction of restorations from pre-milled blanks and in the design and printing of models.
- See also the compatibility table for detailed information on the compatibility of implant scan bodies and model analogues.

Version changes:

2024-02-14 CAD SUPPORT version 1.1.0:

- MIS Original Scanbody added
- Update EV Connection IO FLO
- Update analogue data

2025-04-24 CAD SUPPORT version 1.1.1:

- Functionality for all Exocad systems (EU)
- Compatibility with Amann Girrbach added